

# PREvention™ HLD8

## HIGH LEVEL DISINFECTANT

Manufactured for:

Virox Technologies Inc.  
2770 Coventry Road  
Oakville, ON CANADA L6H 6R1  
1-800-387-7578 | virox.com



### For medical devices and instruments

Active ingredient: Hydrogen Peroxide ..... 2.0%w/w

**INDICATIONS FOR USE:** Prevention™ HLD8 High Level Disinfectant is a reusable high level disinfectant solution for processing heat sensitive medical devices and instruments, for which heat sterilization is not suitable, when used according to the Directions for Use.

Prevention™ HLD8 High Level Disinfectant is intended for use in manual (bucket and tray) systems made from polypropylene, acrylonitrile-butadiene styrene (ABS), polyethylene, polycarbonate plastics, 304 stainless steel and 316 stainless steel. High level disinfection is achieved when all surfaces have been contacted for eight minutes at a temperature of 20°C / 68°F and a minimum concentration of 1.5% hydrogen peroxide.

NOTE: The use of Prevention™ HLD8 High Level Disinfectant solution in automated endoscope reprocessors (AER) must be part of a validated reprocessing procedure. SEE PRECAUTIONARY STATEMENT BELOW FOR USE IN AN AER.

Medical devices reprocessed in Prevention™ HLD8 High Level Disinfectant solution must first be cleaned according to a validated cleaning procedure or standard, such as the ASTM F-1518 "US Standard Practice for Cleaning and Disinfection of Flexible Fiberoptic and Video Endoscopes Used in The Examination of the Hollow Viscera", CAN/CSA-Z314.8-00 (R2005) "Decontamination of reusable medical devices" or other suitable standard.

According to the Directions for Use, Prevention™ HLD8 High Level Disinfectant achieves high level disinfection when used or reused at ≥ 20°C / 68°F for 8 minutes.

**REUSE PERIOD FOR HIGH LEVEL DISINFECTION:** Prevention™ HLD8 High Level Disinfectant solution has demonstrated disinfection efficacy in the presence of organic soil contamination and microbiological burden during reuse.

Prevention™ HLD8 High Level Disinfectant solution may be reused up to a maximum of 21 days, provided the required conditions of hydrogen peroxide concentration, and temperature exist based upon monitoring described in the Directions for Use. DO NOT rely solely on days in use. The hydrogen peroxide concentration of Prevention™ HLD8 High Level Disinfectant solution during its use-life must be verified before each use with the **Prevention™ HLD8 Solution Test Strip**, which will indicate that the Minimum Recommended Concentration (MRC) of 1.5% hydrogen peroxide has been met.

### GENERAL INFORMATION ON SELECTION AND USE

Choose a germicide with the level of antimicrobial activity that is appropriate for the reusable device. Follow the reusable device instruction for use and standard institutional policies. In the absence of complete instructions, use the following process:

For medical devices and instruments, determine whether the reusable device to be reprocessed is a critical or semi-critical device.

Critical Devices are defined as devices that come in contact with sterile tissue, blood stream or recirculating body fluids. These include surgical instrumentation, implantable devices and laparoscopes. Sterilization is required. DO NOT use Prevention™ HLD8 High Level Disinfectant solution as the final treatment prior to use.

Semi-critical devices are devices that come in contact with damaged skin or intact mucus membranes such as the respiratory tract. These include respiratory equipment, surgical mirrors, flexible endoscopes. High Level Disinfection is required.

### MICROBIAL EFFECTIVENESS

The following indicates the spectrum of activity as demonstrated by testing of Prevention™ HLD8 High Level Disinfectant solution using prescribed test methods.

#### VEGETATIVE ORGANISMS:

*Pseudomonas aeruginosa* (ATCC 15442)  
*Staphylococcus aureus* (ATCC 6538)  
*Salmonella choleraesuis* (ATCC 10708)

#### VIRUSES:

Poliovirus, Chat strain type 1 (ATCC VR-1562)  
Herpes Simplex Type 1  
Adenovirus Type 5

#### SPORES:

*Bacillus subtilis* (ATCC 19659)  
*Clostridium sporogenes* (ATCC 7955)

#### MYCOBACTERIA:

*Mycobacterium terrae* (ATCC 15755)

#### FUNGI:

*Trichophyton mentagrophytes* (ATCC 9533)

### MATERIAL COMPATIBILITY

Prevention™ HLD8 High Level Disinfectant solution has been tested and found to be compatible with the materials shown below.

#### METALS:

Mild Steel  
Gold Plated Steel  
Chrome Plated Steel  
302 Stainless Steel  
304 Stainless Steel  
316 Stainless Steel  
410 Stainless Steel

#### PLASTICS:

High Density Polyethylene (HDPE)  
PTFE (Teflon®)  
Polystyrene  
Polystyrene  
Polycarbonate (Lexan®)  
Polypropylene  
Acrylic  
Polyvinyl Chloride (PVC)  
Acrylonitrile Butadiene Styrene (ABS)  
Nylon

#### ELASTOMERS:

Neoprene  
EPDM 42  
Silicone Rubber  
Viton-A  
Polyurethane  
Natural Rubber (Red)

### PRECAUTIONARY STATEMENTS

#### KEEP OUT OF REACH OF CHILDREN.

#### WARNING

**CONTAINS:** Hydrogen Peroxide. Causes mild skin irritation. Causes eye irritation. Wash hands thoroughly after handling. **IF IN EYES:** Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. If eye irritation persists: Get medical attention. If skin irritation persists: Get medical attention.

Avoid contamination of food in the use and storage of this product. Do not store in food processing areas. Store in dry, well ventilated area away from chemicals, direct light, heat or open flame. Do not mix with other cleaning or disinfecting products.

Refer to the SDS for additional information.

Follow Bloodborne Pathogens Universal Precautions, as defined by authorities such as OSHA, WHO, CSA or other appropriate agency, when handling and cleaning soiled devices. When disinfecting devices, use gloves of an appropriate type, and length, and use proper eye protection. Natural or butyl rubber, nitrile or neoprene gloves are recommended.

Contaminated reusable medical devices and instruments, MUST BE THOROUGHLY CLEANED prior to disinfection, since residual contamination with soil or lubricants will decrease the effectiveness of the disinfectant.

The user MUST adhere to the Directions for Use, as modifications to the Directions for Use may affect the safety and effectiveness of the disinfectant.

Use **Prevention™ HLD8 Solution Test Strip** to confirm hydrogen peroxide concentration before each use. Follow the test strip's Directions for Use provided with each box of **Prevention™ HLD8 Solution Test Strips**. The use of Prevention™ HLD8 High Level Disinfectant solution in automated endoscope reprocessors must be part of a validated reprocessing protocol. The validated AER reprocessing protocol must include solution reuse recommendations. The contact conditions for high level disinfection in the AER must be at a minimum of 20°C for 8 minutes. Use the **Prevention™ HLD8 Solution Test Strip** to monitor hydrogen peroxide concentration before each cycle.

The use of Prevention™ HLD8 High Level Disinfectant with semi-critical devices must be part of a validated rinsing procedure as provided by the device manufacturer. See "DIRECTIONS FOR USE Rinsing Instructions" below for important information.

#### NOTICE TO USER:

Please ensure the instruments, equipment or device material to be disinfected is compatible with Prevention™ HLD8 High Level Disinfectant. If you are uncertain as to material composition of your item, confirm with the manufacturer before proceeding.

**CAUTION:** Corrosive to copper, brass, tungsten carbide, Monel S, silver, chromium-plated brass and nickel-plated steel. May be corrosive to aluminum under prolonged immersion.

### DIRECTIONS FOR USE

**CLEANING/DECONTAMINATION:** Blood and other body fluids must be thoroughly cleaned from hard, non-porous surfaces of medical devices and instruments before application of the disinfectant or sterilant. Follow the instrument manufacturer's instructions for disassembly and cleaning. Perform all necessary leak tests as prescribed by the instrument manufacturer prior to immersion of the instrument in the Prevention™ HLD8 High Level Disinfectant solution. Following thorough cleaning and rinsing, rough dry all instrumentation prior to immersion in Prevention™ HLD8 High Level Disinfectant solution.

NOTE: Additional fluid carried by devices into the Prevention™ HLD8 High Level Disinfectant Solution may lead to increased dilution.

**Prevention™ HLD8 High Level Disinfectant PREPARATION:** NO ACTIVATION OR DILUTION IS REQUIRED. Ensure that the Prevention™ HLD8 High Level Disinfectant solution is within expiration. Record the date the original container was opened on the Prevention™ HLD8 High Level Disinfectant container label, or in a log book. After opening, the solution remaining in the original container may be stored for up to 90 days (provided the 90 days does not extend past the expiration date on the container) until used. Always store remaining solution in its original, closed container.

Pour the desired amount of Prevention™ HLD8 High Level Disinfectant solution from its original container into a secondary container (e.g. soak bucket, AER). Label and record the product name, date dispensed, and expiration date of the solution in the secondary container. The solution in the secondary container can be used for a period up to 21 days. The solution must be discarded after 21 days or sooner as dictated by the **Prevention™ HLD8 Solution Test Strip**. The expiration date of the Prevention™ HLD8 High Level Disinfectant solution CANNOT be extended even if the test strip indicates passing results.

NOTE: It is recommended to cover the secondary container to prevent spillage or extraneous contamination of the solution.

**HIGH LEVEL DISINFECTION:** Place pre-cleaned, rinsed and dried instruments in solutions of undiluted Prevention™ HLD8 High Level Disinfectant. Follow the instrument manufacturer's instructions for reprocessing and the instructions below to ensure that all of the instrument's surfaces have been fully contacted with Prevention™ HLD8 High Level Disinfectant solution. Once the instrument has been immersed and all surfaces in contact with the disinfectant solution, soak the instrument for 8 minutes at 20°C. Monitor the time the instrument is in contact with the solution using a timer.

1. Confirm that the **Prevention™ HLD8 High Level Disinfectant** solution meets the minimum recommended concentration using the Prevention™ HLD8 Solution Test Strip. Follow the test strip instructions for use and interpretation. If the indicator strip indicates that the concentration is not acceptable, discard the solution and DO NOT process instrumentation.
2. After devices are pre-cleaned, rinsed and excess moisture is removed, place in the solution of undiluted Prevention™ HLD8 High Level Disinfectant. Ensure that all lumens are filled with fluid. Follow the instrument manufacturer's instructions for flowing the lumens.
3. Set a timer for eight minutes and allow the device to be immersed in the solution for the entire eight minute period of time.
4. Upon completion of eight minutes immersion, remove the instrument from the secondary container. Follow the instrument manufacturer's instructions for draining the lumens prior to rinsing. Follow the rinsing instructions below.

NOTE: Prevention™ HLD8 High Level Disinfectant may not be suitable for all instruments. **SEE NOTICE TO USER.**

#### RINSING INSTRUCTIONS

5. Following removal from Prevention™ HLD8 High Level Disinfectant solution, thoroughly rinse the medical device by immersing it completely in water. Use sterile water or potable water as required by facility policies. Refer to sections "STERILE WATER RINSE" and "POTABLE WATER RINSE" for further recommendations.
6. Keep the instrument or medical device immersed for a minimum of 1 minute in duration, unless longer is specified by the instrument manufacturer.
7. Manually flush all lumens with large volumes of water (at minimum 100 mL) unless otherwise noted by the instrument manufacturer.
8. Remove device and discard rinse water. Do not reuse the water for rinsing or any other purpose.

NOTE: Refer to the instrument manufacturer's instruction for additional rinsing instructions.

**STERILE WATER RINSE:** The following instruments or medical devices should be rinsed with sterile water, using aseptic techniques when rinsing and handling:

- Devices intended for use in normally sterile areas of the body.
- Devices intended for use in known immunocompromised patients based on institutional procedures (e.g., high risk population served).
- When practical, sterilize bronchoscopes, due to risk of contamination from potable water supply. Although microorganisms in this type of water system are not normally pathogenic in patients with healthy immune systems, AIDS patients or other immunocompromised individuals may be placed at high risk of infection by these opportunistic microorganisms.

**POTABLE WATER RINSE:** For all other devices, a sterile water rinse is recommended when practical. Otherwise, potable tap water rinse is acceptable. When using potable water for rinsing, the user should be aware of the increased risk of recontaminating the medical device with microorganisms which may be present in potable water supplies.

#### TREATED WATER SOURCES:

Water treatment systems, such as softeners or de-ionizers, may add microorganisms to the treated water to the extent that microbial content of the water at the point of use could exceed that of the pretreated drinking water. To ensure proper water quality, adherence to maintenance of the water treatment system is recommended.

The use of a bacterial retentive (0.2 micron) filter system or ultraviolet (UV) systems may eliminate or greatly reduce the amount of these waterborne bacteria from the potable water source. Contact the manufacturer of the filter or UV system for instructions on preventive and maintenance and periodic replacement of the filter to avoid colonization or formation of biofilms in the filter.

#### DRYING

A device that is not completely dried provides an ideal situation for rapid colonization of bacteria. As these waterborne bacteria are highly resistant to drying, rapid drying will avoid possible colonization but may not result in a device free from these bacteria. A final rinse using a 70% isopropyl alcohol solution (followed by medical air purge if within lumens) may be used to speed the drying process. Follow facility policies. Consult the instrument manufacturer's instructions for recommended rinsing and drying procedure.

**POST-PROCESSING HANDLING AND STORAGE OF REUSABLE DEVICES:** Disinfected reusable devices are either to be immediately used, or stored in a manner to minimize recontamination. Refer to the reusable device manufacturer's labeling for additional storage and/or handling instructions.

#### EMERGENCY AND TECHNICAL PRODUCT INFORMATION

For further hazard information refer to Safety Data Sheet. Emergency, safety, or technical information about Prevention™ HLD8 High Level Disinfectant solution can be obtained from Virox Technologies at (800) 387-7578, or at www.virox.com.

#### STORAGE CONDITIONS AND EXPIRATION DATE

Prevention™ HLD8 High Level Disinfectant should be stored in its original sealed container at controlled room temperature 15 – 30°C (59 – 86°F). Do not store in food processing areas. Store in dry, well ventilated area away from chemicals, direct light, heat or open flame. Do not mix with other cleaning or disinfecting products. The expiration date of the Prevention™ HLD8 High Level Disinfectant is found on the immediate container.

#### DISPOSAL

Check state and local disposal regulations. Discard residual solution into drain. Flush drain thoroughly with water. Recommended to triple rinse the container with water, and dispose in accordance with federal, state, provincial and municipal requirements. Do not re-use the empty container.

#### HOW SUPPLIED

Catalog Number  
25505

Description  
1.06 Gal (4L) container

Case Contains  
4 containers/case

**Required Consumables: Prevention™ HLD8 Solution Test Strip**